<b>LABiTec -Control Plasma 2- (English)</b>	<b>CE</b>
Abnormal Coagulation Time	

<b>REF</b>	<b>210-13-000-00</b>	<b>CONT</b>	<b>10x1mL</b>
------------	----------------------	-------------	---------------

**INTENDED USE**

The **LABiTec Control Plasma 2** may be used as an abnormal control when measuring PT, APTT and Fibrinogen<sup>(1)</sup>. The **LABiTec Control Plasma 2** is used in the same manner as a fresh plasma pool. The control plasma should be used to gauge internal factors in each laboratory’s system.

**REAGENT**

The **LABiTec Control Plasma 2** is prepared from frozen pool of citrated plasma and is adjusted to give slightly abnormal test results. It is then buffered and lyophilized to ensure stability of all plasma constituents.

**PRECAUTIONS**

- Do not ingest.
- Avoid contact with skin, eyes or clothing.
- WARNING: POTENTIAL BIOHAZARDOUS MATERIAL**  
The source material for this product has been tested and found negative for the presence of HIV and HCV antibodies as well as Hepatitis B Surface Antigen by approved test methods. However, no known test method can offer assurance that products derived from human blood are free of infectious agents. Therefore, handle this material observing the same safety precautions employed when handling any potentially infectious material.

**REAGENT PREPARATION**

- Reconstitute the **LABiTec Control Plasma 2** with exactly 1mL of purified water.
- Replace the stopper and gently mix the vial to thoroughly disperse the contents. Let stand at room temperature for no less than 30 minutes before use to assure complete rehydration of the contents.

**STORAGE AND STABILITY**

The **LABiTec Control Plasma 2** is stable according to the following specification.

Lyophilized	2 – 8°C	Refer to vial label
Reconstituted	18 to 26°C	8 hours
Frozen*	-20°C	2 months
Frozen*	-70°C	6 months

\* Quick thaw frozen controls and test immediately. The controls must not have any contact with glass. Do not incubate controls at 37°C for longer than 5 minutes to avoid the loss of factors V and VII. Loss of factor V and VII can prolong the PT reading.

**PROCEDURE**

The reconstituted **LABiTec Control Plasma 2** is tested in the same manner as freshly drawn citrated patient plasma in Prothrombin time test, activated partial thromboplastin times and Fibrinogens. Refer to the appropriate product inserts for test specific instructions.

**LIMITATIONS**

The results obtained depend on several factors strongly associated with instrument, types of reagent and laboratory to laboratory variations<sup>(2-4)</sup>. Each laboratory should establish an expected range for the particular instrument-reagent system used.

**PERFORMANCE CHARACTERISTICS**

Influences such as reagent type, ISI value of the PT reagent, methodology, instrumentation and technique contribute to variation in test results. Each laboratory should establish its own acceptance ranges with each new lot of plasma control. The **LABiTec Control Plasma 2** will typically yield results within the range specified in the following table, for LABiTec’s PT, APTT and Fibrinogen assays.

<b>Coagulation Test</b>	<b>Abnormal Clotting Time (seconds)</b>
Prothrombin Time (PT) ISI ~ 1.0	16 – 32
Activated Partial Thromboplastin Time (APTT)	40 – 70
Fibrinogen g/L	0.6 – 1.3

The coefficient of variation (CV) for Prothrombin time (PT) and activated partial thromboplastin time (APTT) tests performed on the **LABiTec Control Plasma 2** has been shown to be less than 5 % in intra-laboratory studies. However, precision characteristics will vary depending on the instrumentation and reagent system used.

The results are shown in the following table:

	<b>PT Precision</b>	<b>APTT Precision</b>	<b>Fibrinogen Precision</b>
Within-run (n=20)	± 3.0 <span> </span> % CV	± 1.5 <span> </span> % CV	± 1.5% CV
Day to Day (5 days)	± 1.8 <span> </span> % CV	± 2.8 <span> </span> % CV	± 3.4% CV

**REFERENCES**

- Morse, E.E., et al., Automated Fibrinogen Determination, AJCP, 55:671-676, 1971.
- Kirkwood, T.B.L. et al., Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559- 568, 1977.
- Goldenfarb, M.D., Reproducibility in Coagulation Assays, AJCP, 55:561-564, 1971.
- Palkuti, H.A., Longberry, J.R., A Precision Study of Coagulation Factor Assay Technics, AJCP, 59:231-235, 1973.

**WARRANTY**

This product is warranted to perform in accordance with its labeling and literature. **LABiTec** disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of **LABiTec’s** products for their specific applications. In no event will **LABiTec** be liable for any consequential damages arising out of aforesaid express warranty.

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LABiTec GmbH An der Strusbek 6 22926 Ahrensburg, Germany www.labitec.com