<b>LABiTec -PT-R- (English)</b>	
Liquid stable recombinant thromboplastin reagent	
System reagents for CoaData and CoaLAB analyzer	
<b>REF</b>	<b>CONT</b>
<b>210-08-010-00</b>	<b>10x2mL</b>
<b>210-08-060-00</b>	<b>10x4mL</b>

#### INTERNATIONAL SENSITIVITY INDEX

**Target ISI Value:** ~1.0 (Actual value refer to table of values attached)

#### INTENDED USE

The **LABiTec PT-R** is a liquid stable thromboplastin reagent containing human recombinant tissue factor (RTF)<sup>(1)</sup>. It is used for the quantitative determination in human citrated plasma of Prothrombin Time (PT) and evaluation of extrinsic coagulation factor assays II, V, VII and X.

#### SUMMARY

The **LABiTec PT-R** is a liposomal preparation that contains human RTF relipidated in a phosfolipid blend combined with calcium chloride, buffer and a preservative. **LABiTec PT-R** is especially suitable for monitoring oral anticoagulant therapy (OAT). The **LABiTec PT-R** is particularly appropriate for extrinsic pathway factor assays; this is due to the fact that RTF does not contain any contaminating coagulation factors. The **LABiTec PT-R** reagent is formulated to be insensitive to therapeutic levels of heparin. In the PT test the addition of the tissue thromboplastin to the patient plasma in the presence of calcium ions initiates the activation of the extrinsic pathway.

#### REAGENT

The **LABiTec PT-R** is a liquid stable thromboplastin reagent containing human recombinant tissue factor (RTF), calcium chloride, buffer and sodium azide as a preservative. The lot number of the reagent is shown on the vial's label.

#### PRECAUTIONS

Do not ingest. Avoid contact with skin, eyes or clothing.

#### REAGENT PREPARATION

**LABiTec PT-R is supplied ready to use. Mix reagent by gentle inversion prior to use.**

#### STORAGE AND STABILITY

The **LABiTec PT-R** is stable according to the following specification. DO NOT FREEZE

Unopened	2 – 8°C	Refer to vial label
Open Vial	2 – 8°C	7 Days
On-board CoaLAB 1000	18 - 26°C	5 Days
On-board CoaData Instruments	37°C	8 Hours

#### SPECIMEN COLLECTION AND PREPARATION

Test plasma should be prepared from citrated whole blood **without** heparin, EDTA or oxalate. To obtain the plasma, carefully mix 1 part sodium citrate solution (0.11 mol/L) with 9 parts venous blood, avoiding the formation of foam. Immediately centrifuge for 10 minutes at approximately 3000 rpm (1500 x g), remove the supernatant plasma and keep at + 15 to + 25°C until use in the test. In the US please refer to the NCCLS Document H21-A2, entitled "Collection, Transport and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays".

#### PLASMA STORAGE:

Stability of Sample Plasma		
At ambient temp.	18 to 26°C	8 hours
Frozen*	-20°C	2 months
Frozen*	-70°C	6 months

\* Quick thaw frozen samples and test immediately. The samples must not have any contact with glass. Do not incubate samples at 37°C for longer than 5 minutes to avoid the loss of factors V and VII. Loss of factor V and VII can prolong the PT reading.

#### MATERIALS PROVIDED

#### LABiTec PT-R

#### MATERIALS REQUIRED BUT NOT PROVIDED

- LABiTec Control Plasmas 1 and 2**

#### PROCEDURE

This procedure pertains to manual or semi-automated coagulation systems. Refer to your instrument manual for more detailed instrument specific instructions.

Pre-incubate the <b>LABiTec PT-R</b> to 37°C for at least 10 minutes (water-bath). Maintain the suspension of the reagent by mixing by inversion immediately prior to use.	
Stepwise pipetting and test procedure for manual systems:	
<b>50µl</b>	Sample, control or standard plasma.
<b>60s</b>	Incubation at 37°C (in the measuring block).
<b>100µl</b>	<b>LABiTec PT-R</b> (pre-warmed at 37°C) as start reagent.
The clotting time in seconds will be recorded automatically by the analyzer. Refer to the instrument operator's manual for further information.	

#### QUALITY CONTROL

Reliability of test results should be monitored by use of LABiTec Control Plasma 1-2. For results continuously out of range a local INR calibration<sup>(2,3)</sup> is recommended. A control range should be established by the laboratory to determine the allowable variation in day to day performance for each level of control. Failure to meet quality control specifications should be investigated and resolved.

<b>LABiTec Control Plasma 1</b>	(REF 210-12-000-00)	Normal Values
<b>LABiTec Control Plasma 2</b>	(REF 210-13-000-00)	Abnormal Values

In case the values of Control Plasma are out of the given ranges, check: Calibration, reagent, stability of sample/controls or analyzer.

#### INTERNATIONAL NORMALIZED RATIO (CALCULATION OF RESULTS)

The International Committee for Standardization in Hematology and the International Committee on Thrombosis and Hemostasis have agreed on recommendations for the reporting of Prothrombin Time results as an International Normalized Ratio (INR). The INR is based on the International Sensitivity Index (ISI) of Thromboplastin reagents.

PT reagents are assigned an ISI value by calibration against an International Reference Preparation (IRP 67/40) with an assigned ISI of 1.0. The ISI value for each lot of **LABiTec PT-R** appears on the vial label. Mean Normal PT value appears on each lot specific table of values. The INR is calculated using the following formula:

**INR = (Patient PT / Mean Normal PT)<sup>ISI</sup>**

**ISI** = Lot specific International Sensitivity Index for the Reagent/Instrument system

**Mean Normal PT** = Lot specific mean of the normal range, as determined by each laboratory for the Reagent /Instrument System. It is usually based upon the mean PT plus or minus 2 to 3 standard deviations using 20 or more normal individuals.

#### LIMITATIONS

To prevent spurious results ensure the blood to anticoagulant ratio is 9+1 (1:10). The PT clotting times may be prolonged by substances including corticosteroids, EDTA, oral contraceptives, asparaginase, clofibrate, erythromycin, ethanol, tetracycline and anticoagulants such as heparin and Coumadin. The PT may be shortened by substances including antihistamines, butabarbital, caffeine, oral contraceptives, phenobarbital and vitamin K.

#### EXPECTED VALUES

<b>Analyzer</b>	<b>Method</b>	<b>Normal Values (in seconds)</b>
CoaData	Manual	10 - 15
CoaLAB	Automated	10 - 15

#### TABLE OF VALUES AND CALIBRATION

Each lot of **LABiTec PT-R** is provided with a pre-calibrated table of values for conversion of PT readings to INR and percent PT. The table of values are specifically calibrated for LABiTec's CoaData and CoaLAB analyzers. It contains the **MNPT (Mean Normal PT)** value which represents the 100% value of the reference curve. The recommended calibration should be verified with the **LABiTec Control Plasmas 1 and 2** and it can be used for routine if the measured results are in the provided ranges. Otherwise purchaser is responsible to re-calibrate and determine the suitability of LABiTec's products for their specific applications.

Each laboratory should establish a normal range using individual's representative of its patient population. A new normal range should be established with any change of instrumentation, blood collection techniques, or anticoagulant. A Mean Normal PT range should be reestablished or verified when changing lots of the same reagent. Therapeutic ranges for monitoring oral anticoagulant therapy will vary from laboratory to laboratory. Therefore, it is essential that each laboratory establish relevant PT ranges for its respective patient population.

Abnormal results obtained with a plasma from a patient not on anticoagulant therapy may indicate a factor deficiency or the presence of an inhibitor. The result may also be due to the effects of certain drugs and medications. Additional procedures such as an APTT and mixing studies using factor deficient plasma are usually required.

#### PERFORMANCE CHARACTERISTICS

**Precision:** Within-run precision was assessed using **LABiTec Control Plasmas 1 and 2**, on an automatic and a manual instrument. The results are shown in the following table.

<b>Within-run Precision Results:</b>			
<b>Sample</b>	<b>CoaLAB (Automatic)</b>	<b>CoaData (Manual)</b>	
Level 1	2.06 <span> </span> %		1.93 <span> </span> %
Level 2	3.00 <span> </span> %		3.14 <span> </span> %

**Correlation:** Correlation studies were performed against the PT Recombiplastin reagent of a competitor on the CoaLAB coagulometer. The results are shown in the following table.

<b>Correlation Results</b>			
<b>Regression coefficient</b>	<b>Slope</b>	<b>Intercept</b>	
0.999	1.43	-4.01	

#### REFERENCES

- Tripodi, A etal. (1992) Recombinant tissue factor as substitute for conventional thromboplastin in the Prothrombin Time test. Thrombosis Hemostasis. 67 (1) 42-45.
- Van Rijn, JLM.L. (1989) Correction of Instrument – and Reagent based in Difference in Determination of the International Normalised Ration (INR) for Monitoring Anticoagulant Therapy, Clinical Chemistry. 35 (5) 840-843.
- ICSH/ISTH (1985) Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control. Thrombosis Haemostasis. 53 (1) 155-156
- Lewis, SM. Bain, B.J. Baites, I. (2001) Pratical Haematology Ninth Edition. Churchill Livingstone Publications. 353-354.

#### WARRANTY

This product is warranted to perform in accordance with its labeling and literature. **LABiTec** disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of **LABiTec's** products for their specific applications. In no event will **LABiTec** be liable for any consequential damages arising out of aforesaid express warranty.

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**BIBLIOGRAFIA**

- Tripodi, A etal. (1992) Recombinant tissue factor as substitute for conventional thromboplastin in the Prothrombin Time test. Thrombosis Hemostasis. 67 (1) 42-45.
- Van Rijn, JLM.L. (1989) Correction of Instrument – and Reagent based in Difference in Determination of the International Normalised Ration (INR) for Monitoring Anticoagulant Therapy, Clinical Chemistry. 35 (5) 840-843.
- ICSH/ISTH (1985) Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control. Thrombosis Haemostasis. 53 (1) 155-156
- Lewis, SM. Bain, BJ. Baites, I. (2001) Pratical Haematology Ninth Edition. Churchill Livingstone Publications. 353-354.

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**REFERENCIAS**

- Tripodi, A etal. (1992) Recombinant tissue factor as substitute for conventional thromboplastin in theProthrombin Time test. Thrombosis Hemostasis. 67 (1) 42-45.
- Van Rijn, JLM.L. (1989) Correction of Instrument – and Reagent based in Difference in Determination of the International Normalised Ration (INR) for Monitoring Anticoagulant Therapy, Clinical Chemistry. 35 (5) 840-843.
- ICSH/ISTH (1985) Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control. Thrombosis Haemostasis. 53 (1) 155-156
- Lewis, SM. Bain, BJ. Baites, I. (2001) Pratical Haematology Ninth Edition. Churchill Livingstone Publications. 353-354.

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