LABiTec -Standard Plasma- (English) Calibration Plasma for Fibrinogen and Prothrombin Time (PT) 
--

 210-14-000-00  10x1mL

INTENDED USE

The **LABiTec Standard Plasma** is intended for use in performing calibration curves of Fibrinogen and Prothrombin Time (PT) in %.

The dilutions of Fibrinogen are done with Imidazole Buffer and physiological saline solution is used for PT.

REAGENT

The **LABiTec Standard Plasma** is prepared from a frozen pool of citrated plasma and is buffered and lyophilised to ensure stability of all plasma constituents.

WARNING AND PRECAUTIONS

- Do not ingest.
- Avoid contact with skin, eyes or clothing.
- WARNING: POTENTIAL BIOHAZARDOUS MATERIAL**

The source material for this product has been tested and found negative for the presence of HIV and HCV antibodies as well as Hepatitis B Surface Antigen by approved test methods. However, no known test method can offer assurance that products derived from human blood are free of infectious agents. Therefore, handle this material observing the same safety precautions employed when handling any potentially infectious material.

REAGENT PREPARATION

- Reconstitute the **LABiTec Standard Plasma** with exactly 1mL of purified water.
- Replace the stopper and gently mix the vial to thoroughly disperse the contents. Let stand at room temperature for no less than 30 minutes before use to assure complete rehydration of the contents.

STORAGE AND STABILITY

The **LABiTec Standard Plasma** is stable according to the following specification:

Lyophilized	2 – 8°C	Refer to vial label
At ambient temp.	18 to 26°C	8 hours
Frozen*	-20°C	2 months
Frozen*	-70°C	6 months

* Quick thaw frozen plasma and test immediately. The plasma must not have any contact with glass. Do not incubate plasma at 37°C for longer than 5 minutes to avoid the loss of factors V and VII. Loss of factor V and VII can prolong the PT reading.

PROCEDURE FOR USE

The reconstituted **LABiTec Standard Plasma** is tested in the same manner as freshly drawn citrated patient plasma in prothrombin time, thrombin time, activated partial thromboplastin time tests and fibrinogen assays. Refer to the appropriate product inserts for test specific instructions.

LIMITATIONS

The results obtained depend on several factors strongly associated with instrument, types of reagent and laboratory to laboratory variations ⁽²⁻⁴⁾. Each laboratory should establish an expected range for the particular instrument-reagent system used.

Coagulation Test	Normal Clotting Time (seconds)
Prothrombin Time (PT)	10 – 15
Fibrinogen g/L	1.8 – 4.0

Recommendation to perform a calibration curve:

Dilution	Standard Concentration	Pipette into tubes	add LABiTec Standard Plasma
1 + 4	6.00 g/L	400µl Imidazole Buffer	+ 100µl
1 + 9	3.00 g/L	250µl Imidazole Buffer	+250µl (1 + 4 diluted)
1 + 19	1.50 g/L	250µl Imidazole Buffer	+250µl (1 + 9 diluted)
1 + 39	0.75 g/L	250µl Imidazole Buffer	+250µl (1+19 diluted)

Standard Concentration	Pipette into tubes
100 %	500µl LABiTec Standard Plasma
50 %	250µl saline + 250µl 100% LABiTec Standard Plasma
25 %	250µl saline + 250µl 50% LABiTec Standard Plasma
12.5 %	250µl saline + 250µl 25% LABiTec Standard Plasma

REFERENCES

- Morse, E.E., et al., Automated Fibrinogen Determination, AJCP, 55:671-676, 1971.
- Kirkwood, T.B.L. et al., Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559- 568, 1977.
- Goldenfarb, M.D., Reproducibility in Coagulation Assays, AJCP, 55:561-564, 1971.
- Palkuti, H.A., Longberry, J.R., A Precision Study of Coagulation Factor Assay Technics, AJCP, 59:231-235, 1973.

WARRANTY

This product is warranted to perform in accordance with its labeling and literature. **LABiTec** disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of **LABiTec's** products for their specific applications. In no event will **LABiTec** be liable for any consequential damages arising out of aforesaid express warranty.

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LABITec GmbH