LABiTec -Fibrinogen Kit- (English) <div>For quantitative determination of Fibrinogen</div> <div>System reagents for CoaData and CoaLAB analyzer</div>	
--	--

REF	210-11-000-00	CONT	5x 2ml 4x 25ml 1x 11ml
------------	----------------------	-------------	---

INTENDED USE

The **LABiTec Fibrinogen Kit** is an *in vitro* diagnostic assay intended for quantitative determination of fibrinogen in plasma. The **LABiTec Fibrinogen Kit** is based on the Clauss method ⁽¹⁾. The Clauss method measures the rate of fibrinogen to fibrin conversion in the presence of excess thrombin and has been shown to be rapid, sensitive and precise.

SUMMARY

Thrombin converts soluble fibrinogen into insoluble fibrin, which when cross-linked becomes the fibrin clot as the last step in the coagulation cascade. Fibrinogen is an acute-phase reactant protein in that the concentration rises sharply in response to many different physiological stimuli such as tissue inflammation or injury. High fibrinogen levels are associated with atherosclerotic cardiovascular disease and with the occurrence of myocardial infarction and stroke. Other conditions in which fibrinogen is elevated are cancers of the stomach, breast, or kidney, and inflammatory disorders like rheumatoid arthritis. Reduced fibrinogen levels are prevalent in liver disease, prostate cancer, lung disease, bone marrow lesions, malnourishment, and disseminated intravascular coagulation. Other conditions of deficient fibrinogen are congenital afibrinogenemia, hypofibrinogenemia, and dysfibrinogenemia.

REAGENTS

- The **LABiTec Thrombin Reagent (R1)** is a lyophilized preparation of bovine thrombin, approximately 100 NIH U/mL, buffer, stabilizers and preservative. Reconstitute individual vials with 2 mL of the **Kaolin Suspension (R3)** provided. Allow to stand at room temperature for 30 minutes before use.
- The **LABiTec Imidazole Buffer (R2)** is ready for use. Store at 2-8°C. Use until expiration date printed on bottle. Avoid contamination. Used for sample dilution.
- The **LABiTec Kaolin Suspension (R3)** is ready for use. **Shake well before use!**

STORAGE AND STABILITY

The **LABiTec Fibrinogen Kit** is stable according to the following specification.

Lyophilized	2 – 8°C	Refer to vial label
Reconstituted	2 – 8°C	4 days
On-board CoaLAB 1000	18 - 26°C	8 Hours
CoaData Instruments	18 - 26°C	8 Hours

PRECAUTIONS

Do not ingest. Avoid contact with skin, eyes or clothing. Fibrinogen Plasma Controls whenever used with this assay, are potentially biohazardous materials. Source materials from which plasma controls are manufactured were found negative for HBsAg and for antibodies against HCV, HIV-1 and HIV-2 using approved methods; however, no test method can offer complete assurance that infectious agents are absent. As with all materials of human origin, plasma controls should be handled as a potentially infectious material.

HAZARD STATEMENTS

Imidazole: H360D May damage the unborn child.

SPECIMEN COLLECTION AND PREPARATION

Test plasma should be prepared from citrated whole blood **without** heparin, EDTA or oxalate.

To obtain the plasma, carefully mix 1 part sodium citrate solution (0.11 mol/L) with 9 parts venous blood, avoiding the formation of foam. Immediately centrifuge for 10 minutes at approximately 3000 rpm (1500 x g), remove the supernatant plasma and keep at + 15 to + 25°C until use in the test. In the US please refer to the NCCLS Document H21-A2, entitled "Collection, Transport and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays".

PLASMA STORAGE

Stability of Sample Plasma		
Undiluted at ambient temp.	18 to 26°C	8 hours
1+9 diluted at ambient temp.	18 to 26°C	1 hour
Frozen*	-20°C	2 months
Frozen*	-70°C	6 months

* Quick thaw frozen samples and test immediately. The samples must not have any contact with glass. Do not incubate samples at 37°C for longer than 5 minutes to avoid the loss of factors V and VII. Loss of factor V and VII can prolong the PT reading.

MATERIALS PROVIDED

- LABiTec Thrombin Reagent (R1)** 5x 2mL
- LABiTec Imidazole Buffer (R2)** 4x 25mL
- LABiTec Kaolin Suspension (R3)** 1x 11mL

MATERIALS REQUIRED BUT NOT PROVIDED

- LABiTec Control Plasmas 1 and 2**
- LABiTec Standard Plasma** (REF 210-14-000-00)

PROCEDURE

This procedure pertains to manual or semi-automated coagulation systems. Refer to your instrument manual for more detailed instrument specific instructions.

Reconstitute the LABiTec Thrombin Reagent (R1) with LABiTec Kaolin Suspension (R3) .	
Stepwise pipetting and test procedure for manual systems:	
100µl	1+9 diluted sample, control or standard plasma with (R2) .
60s	Incubation at 37°C (in the measuring block).
50µl	LABiTec Thrombin (R1) (ambient temp.) reconstituted with LABiTec Kaolin Suspension (R3) as start reagent.
The clotting time in seconds will be recorded automatically by the analyzer. Refer to the instrument operator’s manual for further information.	

CONVERSION TABLE

Each lot of **LABiTec Fibrinogen Kit** is provided with a pre-calibrated conversion table for determining fibrinogen concentrations corresponding to FIB clotting times of patient samples. These calibrations are specific to LABiTec’s CoaLAB and CoaData Instruments. Users of other instruments should perform own calibrations as outlined above.

QUALITY CONTROL

Reliability of the calibration curve should be monitored within each run using normal and abnormal fibrinogen control plasmas. Each laboratory should establish a control range to determine the allowable variation in day-to-day performance of each control plasma. The standard curve should be prepared on a monthly basis or when a new lot of reagent is used to assure proper performance.

LABiTec Control Plasma 1	(REF 210-12-000-00)	Normal Values
LABiTec Control Plasma 2	(REF 210-13-000-00)	Abnormal Values

In case the values of Control Plasma are out of the given ranges, check: Calibration, reagent, stability of sample/controls or analyzer.

LIMITATIONS

If the clotting time of the 1:10 dilution of test plasma exceeds the clotting time of the last dilution point on the calibration curve, make a 1:5 dilution of test plasma and repeat the assay. Multiply the resulting value from the curve by 5 instead of 10 to allow for the different dilution factor. This will give the final concentration of the undiluted patient plasma.

If the clotting time of the 1:10 dilution of test plasma is shorter than the clotting time of the last dilution point on the calibration curve, make a 1:20 dilution of test plasma and repeat the assay. Multiply the resulting value from the curve by 20 instead of 10 to allow for the different dilution factor. If other dilutions are tested, the value obtained should be multiplied by the appropriate dilution factor.

The lowest recommended dilution is 1:3. Undiluted plasma cannot be tested because interfering substances and inhibitors may affect the accuracy of the results. Results are not significantly affected by the usual therapeutic levels of heparin up to 3.0 U/mL as found in anticoagulated patients. Prolonged clotting times will result at approximately 5.0 U/mL in the undiluted patient sample. Fibrin degradation products (FDP) may inhibit the thrombin action on fibrinogen and fibrin polymerization. In samples with normal fibrinogen levels, FDP has minimal effect; however, in samples with fibrinogen concentrations below 1.5 g/L and FDP concentrations greater than 100 µg/mL, the assay may be increasingly inhibited. Further dilution of the test plasma will reduce this interference.

EXPECTED VALUES

The normal range for fibrinogen levels in human plasma is considered to be 1.8 - 4.0g/L. Each laboratory should establish its own mean normal and normal range because of variances among different laboratories.

PERFORMANCE CHARACTERISTICS

Precision: Within-run precision was assessed using **LABiTec Control Plasmas 1 and 2**, on an automatic and a manual instrument. The results are shown in the following table.

Within-run Precision Results:		
Sample	CoaLAB (Automatic)	CoaData (Manual)
Control Plasma 1	4.50 %	1.68 %
Control Plasma 2	3.47 %	1.50 %

REFERENCES

- Clauss, A. (1957) Rapid physiological coagulation method in determination of fibrinogen. ActaHaematologica. 17.237
- Lewis, SM. Bain, BJ. Baites, I. (2001) Practical Haematology Ninth Edition. Churchill Livingstone Publications. 355-357
- Guidelines On Fibrinogen Assays. (2003) British Journal of Haematology. 121, 396–404

WARRANTY

This product is warranted to perform in accordance with its labeling and literature. **LABiTec** disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of **LABiTec's** products for their specific applications. In no event will **LABiTec** be liable for any consequential damages arising out of aforesaid express warranty.

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