<p>CoaPT (English)</p> Prothrombin Time	 		
Reagents suitable for CoaData and CoaLAB analyzer			
[REF]	210-24-010-00	[CONT]	8x5mL

INTENDED PURPOSE (see References 1, 4, 7)

The Prothrombin Time (PT) is a key test to investigate the function of coagulation factors belonging to the extrinsic pathway (Factors II,V,VII,X). An elongated PT generally indicates a low level of one or more factors of the extrinsic pathway, due to inherited or acquired bleeding disorders, vitamin K deficiency, liver disease, fibrinolysis, disseminated intravascular coagulation (DIC), oral anticoagulant therapy (OAT). The anticoagulant drugs reduce the activity of vitamin K dependent coagulation factors (factor II, VII, IX, X, Protein C, Protein S) increasing the Prothrombin Time. The ability of the blood to form a fibrin clot via the extrinsic pathway requires thromboplastin, calcium, factors I, II, V, VII and X. So, the Factors involved in the "intrinsic" coagulation pathway are bypassed; it means that the deficiencies of Factors VIII, IX, XII are not detected by any of PT Reagent. For this purpose, it is used the aPTT test.

PRINCIPLE

The **CoaPT** Reagent is made up by a thromboplastin preparation extracted from rabbit brain, containing an optimized quantity of calcium ions which makes the Reagent highly sensitive to the deficiency of "extrinsic" pathway Factors, II, V, VII, X. This makes it particularly effective in monitoring oral anticoagulant therapy (OAT). The performance of this reagent PT are similar to those of the WHO Reference Material (NIBSC RBT/16).

PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However, it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product.
Material safety data sheet should be supplied on request.

ATTENTION!

A) The reagent can be used with manual, mechanical, photometric and nephelometric clot detection systems. The automated determinations must be performed according to specific instructions attached to the instrument used.
B) Very deep attention must be given to interfering substances: certain drugs and other substances may influence levels of PT or PT assay (see References 2).
C) The reagent must be used ONLY for the intended purpose, by expert and trained people and according to good laboratory practice.
D) The clinical diagnosis cannot be done correctly using the result of only one test but have to be done integrating critically the results of different laboratory tests and clinical data.
E) A series of factors, such as ambient temperature, the temperature of the working reagents, the accuracy of the washings, the type of coagulometer, can affect the performance of the test.
F) When changing the reagent lot, it is important to ensure that the new ISI value (see below) is adopted and a new average value of the PT of normal samples is determined.
G) The calibration curve has always to be repeated at each change of reagent lot and/or calibrator.
H) When using reagents, standard laboratory precautions should be observed. All calibrators and controls are human samples and therefore potentially infectious. All appropriate protective measures must be taken to avoid any potential biohazards.

REAGENTS

Components of the kit: **[REF] 210-24-010-00**
R1 - THROMBOPLASTIN L 8 x 5 mL
Thromboplastin from rabbit brain
Excipients and stabilizers
CaCl2 >0,1 g/L
NaN3 < 0,1%

STORAGE AND STABILITY

The **CoaPT** Reagent is stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container, if not exposed to heat sources and/or pressure variations. In case the primary container is damaged, it should be disposed of properly. The **CoaPT** is stable according to the following specification.

Unopened	2 – 8°C	Refer to vial label
Open Vial (in original vial)	2 – 8°C	30 days
On-board CoaLAB 1000	18 - 26°C	7 days
For CoaData Instruments	37°C	4 hours

We suggest to stir reagent continuously during the tests. At the end of the working cycles is suggested to store the reagent to 2-8°C in the original vial for better stability.

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right test performances use following kits (see the relative instructions for use (IFU)):

CoaCAL COAGULATION CALIBRATOR **[REF] 210-24-070-00**
CoaCPN CONTROL PLASMA NORMAL **[REF] 210-24-050-00**
CoaCPA CONTROL PLASMA ABNORMAL **[REF] 210-24-060-00**

Each laboratory should establish its own range of control that determines the acceptable variation of the day-to-day performance for each control plasma. The calibration curve must be repeated for each new lot of reagents and / or calibrator.

PREPARATION OF THE WORKING REAGENT

The **CoaPT** Reagent is supplied ready to use. Before use, the container must be mixed several times by inversion to ensure the homogeneity of the reagent. Leaving the bottle unused for a long time a thromboplastin sediment might form. Leave the reagent at 15-25°C for 30 minutes. Mix gently before use: DO NOT shake. NEVER freeze the **CoaPT** Reagent. Close immediately after use. Reagents must be used correctly to avoid contamination. Improper handling releases us from any responsibility.

MATERIALS REQUIRED BUT NOT PROVIDED

CoaCAL COAGULATION CALIBRATOR **[REF] 210-24-070-00**
CoaCPN CONTROL PLASMA NORMAL **[REF] 210-24-050-00**
CoaCPA CONTROL PLASMA ABNORMAL **[REF] 210-24-060-00**

Standard laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Plastic test tubes for sample dilution.

Coagulation tubes.

Stopwatch or timer.

Water bath 37°C.

Coagulometers.

SAMPLES

Collection of samples in accordance with CLSI (NCCLS) (see References 3) using citrated tubes.

ANALYTICAL PROCEDURE CoaData SERIES (semi-automated)
Refer to your instrument manual or available application sheets provided for more detailed instrument specific instructions.

Stepwise pipette scheme for manual systems:	
50µl	Sample, control or standard plasma
60s	Incubation at 37°C (in the measuring block).
100µl	CoaPT
The clotting time in seconds will be recorded automatically by the analyzer. Refer to the instrument operator’s manual for further information.	

ANALYTICAL PROCEDURE CoaLAB SERIES (automated Instrument)
Refer to CoaLAB instrument user manual.

ANALYTICAL PROCEDURE for other MANUAL METHOD (see References 6)
• Preheat the **R1 - THROMBOPLASTIN L** at 37°C.
• Pipette into tube coagulation as follows:
- Plasma 0,100 mL
• Incubate the plasma for 1-2 minutes at 37°C.
- **R1 - THROMBOPLASTIN L** 0,200 mL
• Start the stopwatch at the same time the addition of the reagent. Determine the time of formation of the blood clot (seconds).

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

QUALITY CONTROL

Reliability of test results should be monitored by use of **LABITec** Control Plasma **CoaCPN** and **CoaCPA** control range should be established by the laboratory to determine the allowable variation in day-to-day performance for each level of control. Failure to meet quality control specifications should be investigated and resolved.

CoaCPN CONTROL PLASMA NORMAL **[REF] 210-24-050-00**
CoaCPA CONTROL PLASMA ABNORMAL **[REF] 210-24-060-00**

In case the values of Control Plasma are out of the given ranges, check: Calibration, reagent, stability of sample/controls or analyser.

PREPARATION OF CALIBRATION CURVE

Prepare reagent according to the instructions above. Using plasma COAGULATION CALIBRATOR (**CoaCAL** **[REF]210-24-070-00**) and a physiological solution, prepare a dilution series as follows:

COAGULATION CALIBRATOR	100% 0,5 mL	50% 0,5 mL	25% 0,5 mL
Physiological solution	---	0,5 mL	1,5 mL

Determine the PT values of each calibrator prepared (at least in duplicate) WITHIN 30 minutes. Plot on a regular graph the experimental data obtained, placing on the X axis the reciprocal dilutions and on Y the clotting times obtained (seconds). Obtain the best straight line between these 3 points; apply the individual times of patients on this curve, getting everyone to express the result as % of activity compared to the normal value.

ATTENTION!

The kit is tested on Coagulometer CoaLAB 1000 and CoaData 4004 (both **LABITec**) and ACL 7000 (I.L.). Applications on automatic / semi-automatic or manual coagulometer from other manufacturers may be totally different from what we experienced.

CALCULATION OF RESULTS

Calculate the mean clotting times of the samples and controls duplicate; differences of duplicates ≤10%. Repeat the test if necessary. Patient results can be expressed using the following units: seconds to form clots, clotting time ratio of patients versus the average of the normal clotting times, % of activity (this expression requires the Calibration Curve) and International Normalized Ratio (INR). It is recommended INR in patients receiving anticoagulation therapy.

INR (International Normalized Ratio) and ISI (International Sensitivity Index) (see References 4,5,6,7,8)
The International Committees ICSH / ICTH have issued recommendations for reporting the results of PT tests as INR. The INR is based on ISI value of the PT Reagent. For each lot of PT Reagent is assigned an ISI value by the calibration against the WHO International Reference Material (NIBSC RBT/16). The value of ISI will be posted on the outside of each Reagent kit.

The INR is calculated with the following formula:

INR = (patient PT / mean normal PT)^{ISI}

ISI = Batch Index specific to the Reagent/Instrument system.

Mean Normal PT (MNPT) = average of the normal range, specific lot, that each laboratory must determine for the own system Reagent / Instrument, testing at least 20 normal samples.

REFERENCE VALUES (see References 8)

The results of **CoaPT** are influenced by the method of determination of the clot, including reagent and instrument, which may be different from laboratory to laboratory. Generally, the PT test with normal plasmas originates clotting times between 8 and 16 seconds. However, each laboratory must establish its own normal range by testing individuals representative of the population. Even the therapeutic ranges for monitoring the OAT will vary from lab to lab: each laboratory will establish its own reference range PT values for its patient population. Abnormal plasma obtained from non-OAT patients may indicate the lack of a Factor or the presence of an inhibitor: in these cases, it is useful the aPTT test and a study with plasmas Factors deficient.

EXPECTED VALUES

Analyzer	Method	Normal Values INR
CoaData	Manual	0,80 – 1,20
CoaLAB	Automated	0,80 – 1,20

ANALYTICAL PERFORMANCES

(validate on CoaLAB 1000, CoaData 4004 and ACL 7000 (I.L.))

The performances of the **CoaPT** Reagent have been tested with a CoaLAB 1000, CoaData 4004 and an ACL 7000 Coagulometer. The data, while representing the characteristics of the product, could be different for each laboratory and for different coagulometers.

Method Limitations: (see Reference 2)

Interferences: (see References 2)

Interference test criterion: recovery ± 10% of initial value. No interference found on samples with:

- heparin up to 0,5 U/mL;
- total bilirubin up to 40 mg/dL;
- haemoglobin up to 600 mg/dL;
- lipemia [Intralipid ®] up to 1000 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples. The results obtained are following:

Sample	CoaLAB		CoaData		ACL 7000	
	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%
Human Plasma 1	10.8±0.16	0.76	13.7±0.26	0.94	8.82 ± 0.12	0.72
Human Plasma 2	17.9±0.28	0.79	17.8±0.28	0.80	14.20 ± 0.22	0.79

Run-to-run Precision: determined for 5 days with 20 replications for each day, for 2 samples. The results obtained are the following:

Sample	CoaLAB		CoaData		ACL 7000	
	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%
Human Plasma 1	10.8±0.17	0.81	13.7±0.27	1.00	8.75 ± 0.15	0.84
Human Plasma 2	17.9±0.38	1.07	17.8±0.39	1.10	13.99 ± 0.32	1.16

Accuracy: a group of 46 plasma has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:
Linear regression: y = 0,8254x + 2,164
Correlation coefficient: r = 0,9948 n = 46

WARRANTY

This product is warranted to perform in accordance with its labeling and literature. **LABITec** disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of **LABITec**'s products for their specific applications. In no event will **LABITec** be liable for any consequential damages arising out of aforesaid express warranty.

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

REFERENCES

- Tietz Textbook of Clinical Chemistry, Ed. by C.A. Burtis, Elsevier Saunders Co., St. Louis-MO (2012).
- Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACCC Press (2000).
- CLSI(NCCLS) H21-A5; Collection, Handling, Transport and Storage for Hemostasis. Quick Guide.
- CLSI(NCCLS) H47-A2; One-Stage PT Test and aPTT Test.
- van Rijn, J.L.M.L. et al., Clin. Chem. 35, 840 (1989).
- WHO Tech. Rep. Series N°889, 64 (1999).
- Hirsh, J. et al., Chest, 102, 312S (1992).
- ICSH/ICTH Recommend. for reporting PT in OAT, Acta Haemat. 72, 405 (1984).

Symbols key / Symbolschlüssel / Interpretazione simboli / Clave de los Símboles	
	CE Mark / CE Symbol / Marchio CE / Marca CE
	In Vitro Diagnostic Medical Device / In Vitro Diagnostikum / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Catalogue Number / Bestellnummer / Numero di catalogo / Número de catálogo
	Refer to instruction for use / Bedienungsanleitung beachten / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Biological Risk / Biologisches Risiko / Rischio biologico / Riesgo biológico
	Manufactured by / Hersteller / Fabricante / Fabricante
	Distributor/Händler/Distributore/Distribuidora
	Temperature limitations / Temperaturbegrenzung / Limiti di temperature / Limite de temperatura
	Contents sufficient to "n" tests / Inhalt ausreichend für „n“ Tests / Contenuto sufficiente per numero test / Contenidos suficientes para "n" pruebas
	Lot Number / Chargennummer / Codice del lotto / Código de lote
	Expiration Date (YYYY/MM) / Verwendbar bis / Utilizare entro / Fecha de caducidad
	Contents / Inhalt / Contenuto / Contenido
	Control / Kontrolle / Controllo / Control
	Calibrator / Kalibrator / Calibratore / Calibrador
	Reconstitution Volume / Rekonstitutionsvolumen / Volume di ricostituzione / Volumen de reconstitución
	Respiratory sensitisation, hazard categories 1, 1A, 1B – Germ cellmutagenicity, hazzard 1A, 1B, 2 – Carcinogenicity, hazard categories 1A, 1B, 2 – Reproductive toxicity, hazard 1A, 1B, 2 – Specific Target Organ Toxicity – Single exposure, hazard categories 1, 2 – Specific Target Organ Toxicity – Repeated exposure, hazard categories 1, 2 Aspiracion hazard, hazard category 1 (Health Hazard)
	Acute toxicity (oral, dermal, inhalation), hazard category 4 - Skin irritation, hazard category 2. Eye irritation, hazard category 2 - Skin sensitisation, hazard categories 1, 1A - Specific Target Organ Toxicity — Single exposure, hazard category 3- Respiratory tract irritation. Narcotic effects (Health Hazard)

Rev. 210-24-010-00_04.24

	BEN – BIOCHEMICAL ENTERPRISE S.r.l. - via Toselli, 4 - 20127 Milano – Italy Tel: +39 022871829/93 – Fax: +39 022890853 – www.bensrl.it – info
---	--